

User Guide: Semi-Automatic Temperature Analysis

This user guide explains how to use the model for semi-automatic temperature analysis in infrared videos. The model and the accompanying files required for semi-automatic temperature analysis can be obtained from the shared Google Drive Autonomic innervation. A few simple steps will guide you through performing the task:

1. Download the entire folder *Semi-Automatic Temperature Analysis* and save it locally on your computer.
2. This folder should contain the following files and folder:
 - `SemiAutomaticModel.py`
 - `hand_landmarker.task`
 - Folder: *Data*, containing subfolders: 1. *Videos*, 2. *Processed videos* and 3. *Output*
 - User Guide
3. Ensure you have installed Python and Anaconda or Miniconda
4. Install Tesseract by following the steps on this website:
<https://tesseract-ocr.github.io/tessdoc/Installation.html>
5. Run the following lines in your terminal or command prompt:
 - `Pip install imageio`
 - `Pip install opencv-python`
 - `Pip install xlswriter`
 - `Pip install pytesseract`
6. Store `tesseract.exe` on your computer
7. Insert the 'Path to `tesseract.exe`' in the Python script
8. Place the infrared videos to be analysed in the subfolder *Data*, 1. *Videos*
9. Run the Python code